ADA-ADP

ES | FR | EN | IT | PT | IT

BALANZA DE PRECISI3N

BALANCE DE PR3CISION

PRECISION SCALE

BILANCIA DI PRECISIONE

BALANÇA DE PRECI33O

PR3ZISIONSWAAGE

V.1
28/05/2020

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EN

1. SAFETY INFORMATION	40
2. INSTALLATION	40
3. BEFORE USING THE SCALE	40
3.1 SCALE ENVIRONMENT	40
3.2 UNPACKING	40
3.3 LIST OF MATERIALS	40
4. SCALE SET-UP	40
CONDITIONS OF USE	41
SCALE ASSEMBLY WITH DISPLAY CASE	41
SCALE ASSEMBLY WITH ROUND WEIGHING PAN	41
SCALE ASSEMBLY WITH SQUARE WEIGHING PAN	41
CONNECTING THE ADAPTER	41
MOUNTING OF EXTERNAL DEVICES	41
PRE-HEATING TIME	41
LEVELLING THE SCALE	41
5. DISPLAY AND KEYBOARD DESCRIPTION	42
6. CALIBRATION	43
7. BASIC WEIGHING	43
8. PIECE COUNTER FUNCTION	44
9. PERCENTAGE FUNCTION	44
10. LIMIT FUNCTION	45
10.1 ZONE SETTING CLOSE TO 0 WITH NO ACOUSTIC WARNING	45
10.2 LIMITS MODE ACTIVATION	46
10.3 DELETE THE LIMIT VALUES	46
11. PARAMETER CONFIGURATION	47
12. TRANSMISSION VELOCITY CONFIGURATION	48
13. ENABLING THE WEIGHING UNITS	49
14. RESET TO FACTORY SETTINGS	50
15. DISPLAY ILLUMINATION CONFIGURATION	50
16. WEIGHT LIMIT CONFIGURATION	50
16.1 WARNING MODE SETTING	51
16.2 ACTIVATE THE LIMIT SETTING	51
16.3 CHECKING THE UPPER AND LOWER LIMIT	51
17. OPTIONAL FUNCTIONS	52
18. COMMUNICATION WITH A COMPUTER	52
19. INTERFACE RS-232	53
19.1 DESCRIPTION OF THE PIN	53
20. MAINTENANCE AND TROUBLESHOOTING ROUTINE	53
20.1 TROUBLESHOOTING	53
20.2 CLEANING	53
20.3 CLEANING THE STAINLESS STEEL SURFACE	53
20.4 SAFETY REVIEW	53
21. TROUBLESHOOTING	54
22. TECHNICAL CONDITIONS	54
23. TECHNICAL PARAMETERS	55
24. WARRANTY	55

1. SAFETY INFORMATION

To avoid unnecessary damage to the scale, please pay close attention to the following tips:

- Do not use the scale in hazardous areas.
- Disconnect the scale from the power supply before connecting it to (or disconnecting it from) any other installation.
- Excessive electromagnetic interference may cause the display to read incorrectly. Once any electromagnetic interference is reduced, the scale should operate normally.
- Avoid spilling any liquid on the surface of the scale. Use a slightly damp cloth to clean the scale

2. INSTALLATION

- Make sure that the local voltage is in line with the voltage indicated on the rating pan.
- Pay special attention when using the RS-232 C data output.
- Check that the cable connected matches the diagram provided for RS-232 communication. If this is not the case, disconnect the cable immediately
- If there is any visible damage to the equipment or power cord, please disconnect the scale from the power supply, place it in a safe place and do not use it until it has been installed safely.
- The scale can only be connected to our accessories or optional parts. We shall not be liable for any modifications made by the user to the scale, including improper use of the installation or the use of a power cable other than that supplied with the scale.
- If the scale is not working properly, please contact your dealer.

3. BEFORE USING THE SCALE

3.1 SCALE ENVIRONMENT

This scale requires an environment free of high and low temperatures, corrosion, vibrations, draughts and collision hazards.

3.2 UNPACKING

- Unpack the scale and check that it has not been damaged during transport.
- If there is any external damage, please consult the maintenance and troubleshooting section.
- Save the original packaging of the scale for possible future transport. When packing the scale, remove all cables to avoid unnecessary damage

3.3 LIST OF MATERIALS

- Scale.
- Weighing pan.
- Pan holder.
- Network adapter.
- Instruction manual.


4. SCALE SET-UP

Set up the scale in a suitable place, and avoid the following situations:

- High temperatures and direct sunlight.
- Air currents and vibrations
- Excessive humidity level.


4.1 CONDITIONS OF USE

Do not leave the scale in a very humid place for a long time. If the scale is transported from a place with low temperatures to a place with high temperatures, please leave the scale in the new location for 2 hours before connecting it to the power supply.


SCALE ASSEMBLY WITH DISPLAY CASE

Assemble the parts in the order indicated.


SCALE ASSEMBLY WITH THE ROUND WEIGHING PAN

- Install the pan holder.
- Position the weighing pan.

SCALE ASSEMBLY WITH THE SQUARE WEIGHING PAN

- Install the pan holder.
- Position the weighing pan.

CONNECT THE ADAPTER

- Connect the adapter to the scale.
- Connect the adapter to the scale.


MOUNTING EXTERNAL DEVICES

Before installing any external devices (such as printers or computer connections) please turn off the scale or disconnect the power supply.

PRE-HEATING TIME

To ensure accurate readings, the scale must be pre-heated for 30 minutes before use.

LEVELLING THE SCALE

- Level the scale using the two front adjustable feet until the level bubble is in the centre of the circle.
- Each time the scale is installed in a different location, it must be levelled before use.

5. SCREEN AND KEYBOARD DESCRIPTION

1. Capacity

2. Weighing mode symbols
Count - Parts counter
% - Percentage weight

3. Resolution/Division

4. Weighing units

5. Bubble level

6. Upper / lower limit warning / OK

7. OK indicator (stable reading) in normal mode

OK -	Stable reading.	PCS -	Reading in piece counting mode.
G -	Reading in grammes.	-----	The scale is developing a stable reading.
OZ -	Reading in ounces (1g = 0.03527396200 oz)	UNABLE-	Operation error.
CT -	Carat reading (1g = 5.00000000 ct).	HHHHH-	The object on the weighing pan exceeds the capacity of the scale.
DWT -	Reading in pennies.	LLLLL-	The weighing pan is not installed correctly or has been moved
% -	Percentage function.		

	Calibration / adjustment key		Weighing unit conversion key
	Print/escape key		Tare function key
	Part count function key / confirmation		Menu key
	Percentage function key		On / off key

6. CALIBRATION

PROCEDURE

- 1- Level the scale by adjusting the bubble level and preheat the scale for 30 minutes.
- 2- Press the key  if does not show zero once operation has started.
- 3- Press the  key, and the display will show the capacity of the scale, e.g. 3000.
- 4- Press the  key repeatedly, and the display will show the calibration point within the range, e.g. time: 2000, 1000...
- 5- Choose a calibration value and place the corresponding calibration weight on the weighing pan. Press the  key, and the display will show the message **ACAL**.
- 6- When the calibration is completed, the display will show the weight of the calibration weight.

7. BASIC WEIGHING

PROCEDURE

- 1- Turn on the scale by pressing the  key, and the scale will automatically perform an initialisation process.
- 2- Place the container on the weighing pan.
- 3- Press the  key.
- 4- Place the product inside the container.
- 5- Print the result by pressing the  key.

8. PIECE COUNTER FUNCTION

Symbol on the display: **PCS**

PROCEDURE

- 1- Place the container on the weighing pan.
- 2- Press the  key
- 3- Select the number of pieces to be used for reference, press the  pieces key, and the display will show **"qTy10"**; next, press the  key to display the available numbers of pieces (10, 25, 50, 100). The greater the number of pieces, the more accurate the results.
- 4- Place the selected number of pieces on the pan or inside the container.
- 5- Press the  key, and the display will show the number of pieces used for reference, e.g.: 10
- 6- Load the pieces that remain to be counted, and the display will show the total number of pieces, e.g.: 258 pieces.
- 7- To print the results, press  key
- 8- To return to the normal weighing mode, press  key, and the PCS symbol will disappear.

9. PERCENTAGE FUNCTION

Symbol on the display: **%**

PROCEDURE

- 1- Press the  key
 - 2- Place the reference weight on the pan. When the reading is stable, press  and the display will show **"100,000"** or **"100.00"**, which refers to the 100% accuracy of the scale. The % indicator will appear on the display.
 - 3- Remove the reference weight, and the display will show **"0.000"** or **"0.00"**.
 - 4- Place the object to be weighed on the weighing pan.
 - 5- Wait until the **OK** indicator appears on the display; the display reading indicates the deviation from the reference weight as a percentage.
 - 6- Press the  key, and the scale will return to the normal weighing mode.
- To obtain the percentage of a product in a container, perform the following procedure:
1. Place the empty container on the weighing pan and press  key
 2. Pour the reference weight into the container.
 3. Press the  key once the reading is stable, and the display will show **"100,000"** or **"100.00"**.
 4. Remove the container with the reference weight. Place another container on the pan and press the  key to subtract its weight, then wait until the display reading is stable, at which point the display will show **"0.000"** or **"0.00"**.
 5. Pour the product into the container; the displayed value is the difference in percentage between the product and the reference weight.
 6. Press the  key to return to the normal weighing mode.
 7. If desired, the result can be sent to a printer.

10. LIMITS FUNCTION

To check that the weight of the object is within the predefined range of limits.

Symbols on the display: **LOW / HIGH** (with acoustic warning or OK indicator)

1. Enter the menu mode by pressing the  key
2. Press the  key repeatedly until the message "InSPCT" is displayed.
3. Press the  key
4. To set the upper limit (SET-HI), press the  key twice.
5. Set the decimal point position by pressing the  key repeatedly until the decimal point is in the correct position.
6. Confirm the value with the  key
7. Set the upper limit value as follows:
 1. Press  to increase the value.
 2. Press  to decrease the value.
 3. Press  and the value will blink.
 4. Press  to confirm.
8. Press the  key to set the lower limit value.
9. Select the position of the decimal point by pressing the  key twice.
10. Press the  key repeatedly until the decimal point is in the correct position.
12. Confirm by pressing the  key
13. Set the lower limit value as follows:
 1. Press  to increase the value.
 2. Press  to decrease the value.
 3. Press  and the value will blink.
 4. Press  to confirm.

10.1 ZONE SETTING CLOSE TO 0 WITH NO ACOUSTIC WARNING

Within the **INSPECT mode**:

1. Press the  key until "NoNres" is displayed. Press the  key
3. Set the value of the warning:
 1. Press  to increase the value.
 2. Press  to decrease the value.
 3. Press  and the value will blink.
 4. Press  to confirm.

10.2 LIMITS MODE ACTIVATION.

Within the **INSPECT mode**(once the limits have been set):

1. Press the  key repeatedly until **"ENABLE"** is displayed.
2. Press the  key
3. Press the  key
4. Place the weight on the weighing pan
5. If the display shows **"LOW"** it means that the weight of the object is lower than the value of the lower limit, if it shows **"HI"** accompanied by the acoustic warning it means that the weight of the object is higher than the value of the upper limit, and if the display shows **"OK"** it means that the weight of the object is between the lower and upper limit.
6. If desired, the results can be printed out by pressing the  key
7. Remove the weight from the weighing pan

TO EXIT THE LIMITS MODE

1. Press the  key
2. Press the  key repeatedly until **"INSPECT" IS DISPLAYED**
3. Press the  key
4. Press the  key repeatedly until **"DISABLE"** is displayed.
5. Press the  key, and the scale will exit the weight limit function.

10.3 DELETE THE LIMIT VALUES



1. Press the  key
2. Press the  key repeatedly until **"INSPECT" IS DISPLAYED**
3. Press the  key
4. Press the  key repeatedly until **"CLEAR"** is displayed.
5. Press the  key to delete the entered data.

11. PARAMETER SETTING

To enter the configuration mode, press the  key.


1. Press the  key to select the desired parameter and press the  key to confirm.

2. If you want to exit the configuration mode at any time, press the  key repeatedly until the display shows "ESC", press the  key to confirm.


VEr → FIRMWARE VERSION

INITIATE

Once the weight is stable, the backlighting remains on for the given duration.

ESC leave

12. CONFIGURATION PRINT SPEED

There are three print modes:

STABLE Print: print when the display reading is stable.

INSTANT Print: Print immediately after pressing the  key.

INTERVAL Print: printing at predetermined time intervals.

Note: The print function is separate from the line-feed setting. First you set the print function and then the line-feed number.

To set the print mode, follow the procedure below:

1. Press the  key, then without releasing it press the  key: the display will show **"Print"**
2. Press the  key and the display will show **"STABLE"**, which indicates that the programmed print mode is Print when the display reading is stable.
3. Press the  key to select the desired print mode and confirm by pressing .
4. Press the  key until "InTEr" is displayed, then press the  key and the display will show **"SEC"**.
5. Press the  key repeatedly to display the default interval, then press the  key to confirm. **"0 SEC"** means continuous printing.
6. Press the  key to return to the weighing mode.
7. Once the setting has been made, press the  key repeatedly to display the default line feed. (1-18) Select the appropriate line feed parameter, then press the  key to return to normal weighing mode.

TRANSMISSION VELOCITY CONFIGURATION

1. Press the  key
2. Press the  key and the display will show the message **"bAud"**.
3. Press the  key and the display will show **"300"**.
4. Press the  key repeatedly and the display will show the available transmission speeds. Select the desired speed and confirm with the  key. The display will show **"ParTy"** and the parity check will start.
5. Press the  key and the display will show **"nOnE"** (no check) initially.
6. Press the  key and the display will show the other test modes; **"Odd"** means "odd test" and **"EvEn "** means "always test".
7. Select the desired mode and press  to confirm, and the scale will return to normal weighing mode.

13. ENABLE THE WEIGHING UNITS

The weighing units can be activated or deactivated by the user:

1. Press the  key and the display will show **"PrInT"**.
2. Press the  key repeatedly until the message **"unIT"** is displayed.
3. Press the  key and the display will show **"g yes"** which means that unit **"g"** is activated; press the  key to deactivate unit **"g"** press the  key and the display will show **"g no"**, then press the  key to confirm.
4. Follow the same steps to activate or deactivate the other Oz, Ct or dwt weighing units.

14. RESTORE FACTORY SETTINGS

1. Press the  key and the display will show **"PrinT"**.
2. Press the  key repeatedly until you can see the message **"InITIA"**.
3. Press the  key to confirm; the display will show **"BUSY"** and the scale will return to normal weighing mode.

15. DISPLAY ILLUMINATION SETTING

1. Press the  key and the display will show **"PrinT"**.
2. Press the  key repeatedly until **"bLgHT"** is displayed, then press  to confirm and the display will show **"1 min"** which indicates that the display illumination will go off in 1 minute.
3. Press the  key repeatedly to select the desired time: 1, 2, 3, 5, 10, 15, 30 or 60 minutes.
4. Press the  key repeatedly to select the desired time: 1, 2, 3, 5, 10, 15, 30 or 60 minutes

16. WEIGHT LIMIT CONFIGURATION

Configuration of the upper and lower limit values:

1. Press the  key.
2. Press the  key repeatedly until the message **"InSpCT"** is displayed.
3. Press the  key to confirm, the scale will enter the configuration mode and the display will show **"SET HI"** (upper limit configuration).
4. Press the  key, the display will show **"50"** (value initialization).
5. Press the  key, the display will show **"SET Dp"** (setting of the decimal point position). Press the  key repeatedly to select the decimal position.
6. Press the  key to confirm the decimal position. The display will show the initialisation of the value. Press the  key to increase the value and the  key to decrease it.
7. Enter the desired value, then press  and the entered value will flash. To continue the modification, press the  key; to confirm the value press  and the display will show **"SET HI"**.
8. Press the  key and the display will show the message **"SET LO"** (setting the lower limit value); follow the same steps as described above.

16.1 WARNING MODE SETTING

The scale will display a warning when there is no object on the weighing pan or when the weight of the product is less than the pre-defined lower limit value.

1. Press the  key.
2. Press the  key repeatedly until the message **"InSpCT"** is displayed.
3. Press the  key and the display will show **"SET HI"**.
4. Press the  key repeatedly until **"nOnrES"** is displayed.
5. Press the  key and the display will show **"50"** (the default value warns when the weight on the pan is less than 50% of the lower limit value).
6. Press the  key to increase the value or the  key to decrease it.
7. Enter the desired value, press  and the value will flash.
8. Press the  key to confirm, and the display will show **"SET HI"**.

16.2 ACTIVATE THE LIMIT SETTING

1. Press the  key repeatedly until **"EnABLE"** is displayed.
2. Press the  key to confirm.
3. To deactivate the limit setting, press the  key until **"dISAbL"** is displayed.
4. Press the  key to confirm, and the scale will return to the normal weighing mode.

16.3 CHECKING THE UPPER AND LOWER LIMIT

1. Press the  key and the display will show **"Print"**.
2. Press the  key repeatedly until the message **"InSpCT"** is displayed.
3. Press the  key and the display will show **"SET HI"**.
4. Press the  key and the display will show the value of the upper limit.
5. Press the  key and the display will show **"SET HI"** again.
6. Press the  key and the display will show **"SET LO"**.
7. Press the  key and the display will show the value of the lower limit.

NOTE:

1. To exit the configuration menu at any time, press the  key until the message **"ESC"** is displayed, then press the  key to confirm.
2. To delete the values of the upper and lower limits, press  until the message **"CLEAR"** is displayed, then press the key  and the values of both limits will be set to zero.
3. After making any changes to the upper or lower limit, the weight limit mode must be accessed again.
4. The limit values are stored in the scale. There is no need to reconfigure them.

17. OPTIONAL FUNCTIONS

Hydrostatic weighing hook

On the bottom of the scale there is a hole for installing the weighing hook (optional). (The use of the lower weighing hook is not permitted for balances within Legal Metrology).

1. Open the lower cover of the scale (please turn the scale on its side to open the lower cover, not the other way around).
2. Screw the hook to the right in the lower hole, stop screwing the hook if you feel any resistance.
3. Hanging the object from the hook on a suspension line, for example.
4. If necessary, install a hydrostatic bench to prevent drafts.

18. COMMUNICATION WITH A COMPUTER

The scale's key functions can be used via RS-232. The following commands are available:

Weighing unit conversion.

Percentage function.

Tare function.

Instant print.

Calibration with an external calibration weight.

Part counting function.

Print function.

When the scale is connected to a computer, we recommend using the print function immediately. In response to this command, the scale will transmit any number or message that appears on the display.

The output format is as follows:

A B C D E F G H I J K L M

A - +/- Sign field; does not usually appear when it is a positive number, appears when the reading is negative.

B - G: Number and decimal field; spaces are used when there are fewer than six digits.

H - I: Field of spaces.


J: Unit field; describes the weighing unit of the transmitted weights. The scale will transmit G for grams, O for ounces and C for carats.

K: Stability character; corresponds to the OK symbol on the scale display. S means the reading is stable, a space means the reading is unstable.

L: Return character.

M: Line feed character; indicates the line feed.

19. INTERFACE RS-232

This scale adopts the transmission and reception lines of the RS-232 standard.

The data format is:

- 1 start bit.
- 8 data bits include parity.
- 1 stop bit

19.1 DESCRIPTION OF THE PIN

2 - TXD - - - the scale transmits data.

3 - RXD - - - the scale receives data.

5 - GRD - - - ground signal.

Note: CTS signal modes are not used. The peripheral must have a minimum buffer (15 characters).

We recommend that the cable length be 15 metres and that the load impedance of the connected equipment be between 3000 and 7000 ohms with no more than 2500 pf capacity in parallel.

20. ROUTINE MAINTENANCE AND TROUBLESHOOTING

20.1 TROUBLESHOOTING

Only qualified professionals can carry out repairs to the scale. Serious risks may be incurred if the scale is repaired by unqualified personnel.

20.2 CLEANING

- Turn off the scale and disconnect all cables.
- Avoid spilling liquids into the scale.
- Do not use corrosive cleaners (solvents) under any circumstances.
- Clean the scale with a soft cloth.
- Before cleaning the scale, remove the pan, pan support and ring to prevent damage to the weighing system.

20.3 CLEANING THE STAINLESS STEEL SURFACE

All stainless steel components need to be cleaned regularly. Remove the weighing pan and clean it with a soft cloth or sponge. We recommend using a special cleaner for stainless steel. Clean the stainless steel surface of the scale and

then the weighing pan. Make sure there is no dirt before drying the scale with air. If necessary, smear the stainless steel surface with a special oil as additional protection.

20.4 SAFETY REVIEW

If the scale is not working properly:

- Disconnect it from the network immediately, and keep it out of use.
- Keep it in a safe place and do not use it for the time being.
- Contact the Technical Service or your Distributor for repair.

21. TROUBLESHOOTING

DISPLAY	REASON	SOLUTION
-----	The scale cannot obtain a stable reading or the load cell is faulty.	Contact your Distributor.
HHHHH	The actual weight is 5% higher than the maximum capacity of the scale or the load cell is broken.	Remove the weight from the scale or contact the Technical Service.
LLLLLLL	a) The weighing pan is not installed correctly. b) The load cell is defective	a) Install the pan correctly and press the  key. b) Contact the Technical Service.
NOCAL	Calibration is not responding.	See the calibration section and check that the calibration weight is appropriate.
UNABLE	Missing or incorrect data at the time of switching on the scale.	Read the user's manual.
UNSTABLE READING	The scale environment is unstable. Excessive vibration or draughts, or incorrect pan installation	Place the scale in a more appropriate place. Check and install the pan properly.
NO DISPLAY	There is no power.	Check the scale's current transformer.
Apparent wrong weighing display	The value of the weight differs from that of the weight placed on the weighing pan.	Check that there is no standard tare value

22.- TECHNICAL CONDITIONS

Power, voltage, output frequency	AC-DC Adapter, 220V/110V input, 7.5V output (+15% to -20%) 48-60Hz
Consumption	Maximum = 16W, Average = 8W
Operating temperature range	+ 10°C to +30°C (50° F to 86° F)
Permitted temperature range	+ 5°C to +40°C (41° F to 104° F)
The scale can ensure normal work	+ 5°C to +40°C (41° F to 104° F)

23. TECHNICAL PARAMETERS

DISPLAY	ADA-220	ADA-420	ADA-620	ADP-2200	ADP4200 ADP-6200
Capacity	320g	420g	620g	3200g	5000g
Resolution	0.001g		0.001g	0.01g	0.01g
Tare Range	320g	420g	620g	3200g	5000g
Tare Range	0.001g				0.01g
Response time				2.5s	
Value of the weight for calibration	100 or 200g		100, 200 or 500g	1.2 0 3kg	1.2 o 5kg
Minimum class in weight for calibration	F1 Class		F1 Class	F1 Class	F1 Class
Net weight	4.0kg		2.3kg		
Size of the pan	115mm diameter			160mm diameter	
Lower height of display case				230mm	
Scale dimensions	230x310x330mm			230x310x90mm	

24- GUARANTEE

This scale is guaranteed against all manufacturing and material defects for a period of 1 year from the date of delivery.

During this period, GIROPES S.L. will handle the repair of the scale.

This guarantee does not include damages caused by improper use or overloading of the item.

The warranty does not cover the shipping costs (freight) required to repair the scale.

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marca propiedad de | est une marque de | trade mark propriety of: